

# ***FLOW SELECTOR***



**The Flow Selector allows you more treatment options from one gas source.**

***It's Versatile ...***

The Flow Selector allows you to direct the outlet flow from any standard flowmeter to any *ONE* of the three outlets it provides. Select one of the tubing connections or the male DISS outlet.

***It Saves Time ...***

You can deliver a nebulizer treatment without disconnecting the humidifier.

***It Eliminates Clutter ...***

With the Flow Selector's two (2) optional ports,

you can leave the cannula, humidifier and nebulizer hooked up to the Flow selector for easy treatment whenever needed. No more loose tubing lying around waiting to be used for treatments.

***It's Easy to Use ...***

Simply push the button on the bottom of the unit and rotate the desired connection to the large, "ON" arrow for treatment. When the treatment is completed, push the button again and rotate the primary connection back under the "ON" arrow.


***It's Durable ...***

The Flow Selector has a metal DISS female flowmeter connection and a tough nylon body.

***It's Lightweight and Versatile ...***

Weighing less than 4 oz. it adds little weight to the flowmeter and it operates with any flowmeter that has a DISS male outlet connection.





## FLOW SELECTOR FEATURES:

- Give nebulizer treatments without disconnecting the humidifier.
- Reduces changeover time to seconds.
- No more loose tubing laying around.
- Eliminates tubing nipple adaptors (it has 2 of its own).
- Fits all flowmeters that have DISS male outlet fittings.
- Lightweight.



## ORDERING INFORMATION

Part No. ....PM1000



Precision Medical, Inc.  
 300 Held Drive • Northampton, PA 18067  
 Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285  
 Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240  
 Visit us on the Web: [www.precisionmedical.com](http://www.precisionmedical.com)