

DIGITAL

VAC CHECK



PM760

Accurate... Long Lasting... Compact

State of the Art Design...


The Digital Vac Check gauge incorporates digital microprocessor technology - solid state design with no moving parts.

Compact & Versatile...

The unit is small enough to carry in pocket for easy vacuum level accuracy checks wherever required.

Durable...

Made of rigid ABS shock resistant plastic case with aluminum vacuum hose barb port connection.



SPECIFICATIONS

Dimensions:	Width: 1.63, Height: 2.13, Depth: 1.00
Vacuum Range:	0-760 mmHg
Accuracy:	± 1% over full scale range
Connector:	1/4 in. hose barb
Battery:	Lithium Ion (8 to 10 year life)
Weight:	.15 lb.

ORDERING INFORMATION

Part Number	Description
PM760	Digital Vac Check unit



Precision Medical, Inc.
300 Held Drive • Northampton, PA 18067
Phone: (+001) 610-262-6090 Toll-Free Phone: 1-800-272-7285
Fax: (+001) 610-262-6080 Toll-Free Fax: 1-800-352-1240
www.precisionmedical.com