



LOCKING DEVICE FOR ISOLATION VA


PADLOCK



LOCKED IN
OPEN POSITION



LOCKING DEVICE FOR ISOLATION VALVE



NOTE: PADLOCK NOT SUPPLIED
 ($\frac{9}{32}$ " Padlock is Recommended
 per Device to Operate Properly)

NOTE:
 Locking devices can be ordered as separate items
 or as complete assembly on ball-valves.

Model numbers for locking devices:

1/4 - 1/2"	=V - X - LOCK - N - 0205
3/4 - 1"	=V - X - LOCK - N - 0710
1 - 1/4 - 2"	=V - X - LOCK - N - 1220
2 - 1/2"	=V - X - LOCK - A - 25
3"	=V - X - LOCK - N - 30
4"	=V - X - LOCK - N - 40

Model numbers for Isolation Valve with locking devices:

1/2"	=05
3/4"	=07
1"	=10
1 - 1/4"	=12
1 - 1/2"	=15
2"	=20
2 - 1/2"	=25
3"	=30
4"	=40

85 Fulton Way
 Richmond Hill, Ontario
 L4B 2N4, CANADA

Tel: (905) 764-0800
 Fax: (905) 764-0862

www.amico.com