

## PURELINE M6000 Anesthesia Machine With Oxygen Concentrator

### SAVINGS FOR YOU, SAFETY FOR YOUR PATIENTS

Stop paying for tanks, refills, and deliveries.  
Stop wasting time on tank billing, scheduling,  
management, and maintenance.

Note: Vaporizer not included.  
Sold Separately

### ONE OF A KIND ANESTHESIA MACHINE

The Pureline M6000 anesthesia machine with our integrated oxygen concentrator is one of a kind!

Supera has perfected the technology, added safety features like the built-in E-tank manifold that allows mobile use of the machine without power. It also protects your patient in case of a power failure.

The Pureline oxygen concentrator is UL tested and backed by the industry's best warranty.



**CALCULATE SAVINGS**  
Visit [www.superavet.com/ROI](http://www.superavet.com/ROI) to  
calculate YOUR actual savings!



### Key Features:

- Extracts medical-grade oxygen (USP93%)
- Delivers on demand up to 5 lpm
- Built specifically for anesthesia use
- Quiet operation at 40 dBA
- Energy efficient
- Accepts industry standard vaporizers
- FDA approved, UL tested
- Meets 13485:2016 ISO & 60601-1 IE Standards
- Oxygen flow & purity safety alarms
- Thousands in service since 2008
- Made in the USA

# RETURN ON INVESTMENT IN AS LITTLE AS 12 MONTHS!

By eliminating the cost of oxygen tank rental, purchase, refills, delivery charges and tank management time, the Pureline M6000 will pay for itself quickly.

## BENEFITS OF THE PURELINE M6000

- No tanks, hoses, regulators to buy, store or replace...ever!
- Oxygen is “on demand”—never run out!
- Accepts industry standard vaporizers
- 0-4 lpm precision flow meter
- Easy on/off; controls are large and easily accessible
- Stable “A” frame chassis with a low center of gravity
- Small 18” square footprint saves valuable space
- Includes anesthesia accessories: adult & pediatric unilimb rebreathing circuits, non-rebreathing circuit and heavy-duty Supera rebreathing bags
- Quiet 40 dBA operation (best in the industry)
- Includes one high pressure “E” tank for backup and oxygen flush use (the back up e-tank is not filled by machine)



### CALCULATE SAVINGS

Visit [www.superavet.com/ROI](http://www.superavet.com/ROI) to calculate YOUR actual savings!

Heavy-duty shelf holds all monitor types safely and securely.



Stability by design weight is centered within the “A” frame.

# SUPERA

ANESTHESIA  
INNOVATIONS