

# Start Saving Today With The Pureline® M8000 Anesthesia Machine With Oxygen Concentrator

**Pureline**<sup>®</sup>  
OXYGEN CONCENTRATORS

**M8000**



Stop paying for tanks, refills, and deliveries.

Stop wasting time on tank billing, scheduling and management.

*Instead, get your oxygen for free, the natural way ... from the air!*



## Savings for you, safety for your patients.

The **Pureline M8000** anesthesia machine with our integrated oxygen concentrator is one of a kind!

Supera has perfected the technology, added safety features like the built-in E-tank manifold that allows mobile use of the machine without power. It also protects your patient in case of a power failure.

The Pureline oxygen concentrator is FDA approved (USP 93%), UL tested and backed by the industry's best warranty.

## Key Features

- Extracts medical-grade O2 from room air
- Outputs 8 liters per minute
- Energy efficient
- Built specifically for anesthesia machine use
- Oxygen flow & purity safety alarms
- Safety backup O2 supply
- Accepts industry standard vaporizers
- Virtually maintenance free
- Proudly made in Oregon

P/N M8000 Includes safety pop-off valve

P/N M8000NS Includes regular pop-off valve (not shown)

**Note: Vaporizer not included.  
Sold Separately**

## The Pureline® M8000 Pays For Itself!

By eliminating the cost of oxygen tank rental, purchase, refills, delivery charges and tank management time, the Pureline M8000 will pay for itself very quickly.



**Pureline®**  
OXYGEN CONCENTRATORS

## Return On Investment

can be in as little as 12 months!

Visit [www.superavet.com/ROI](http://www.superavet.com/ROI)

to calculate YOUR actual savings!

### The Benefits of the Pureline M8000

- No tanks, hoses, regulators to buy, store or replace...ever!
- Oxygen is “on demand”—never run out!
- Accepts industry standard vaporizers
- 0-4 lpm precision flow meter
- Easy on/off; controls are large and easily accessible
- Stable “A” frame chassis with a low center of gravity
- Small 18” X 22” footprint saves valuable space
- Includes anesthesia accessories: adult and pediatric unilimb rebreathing circuits, non-rebreathing circuit and heavy-duty rebreathing bags
- Quiet 47 dBA operation (best in the industry)
- Includes one high pressure “E” tank for backup and oxygen flush (the back up e-tank is not filled by machine)



Heavy-duty shelf holds all monitor types safely and securely.



Stability by design  
weight is centered  
within the “A” frame.

Shown with optional vaporizer

**SUPERA** | ANESTHESIA INNOVATIONS

**BROWARD A&C**  
MEDICAL GAS SPECIALISTS

Anesthesia Innovations Since 1993