

Model M1000

Non-Rebreathing “A” Frame Stand

Veterinary Anesthesia Machine

Versatile non-rebreathing machine with built in monitor shelf

Superior stability with our new “A” Frame design stand



Note: Vaporizer & E-tank manifold are not included. Sold Separately

M1000 features:

- Built in stainless steel I.V. pole
- 17 1/2" wide X 9" deep heavy duty monitor shelf
- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- Dual O₂ supply input with check valves
- 15mm female common gas outlet for circuits
- Syringe holder (for inflating E/T cuffs)
- Includes non-rebreathing circuit
- Mounting option for universal control arm for “Bain” circuits (see MA2013 below)



I.V. Pole included

Optional Universal Control Arm

Model: MA2013

