

# SUPERA

ANESTHESIA  
INNOVATIONS

## Model M2000 Deluxe Rebreathing Floor Stand Veterinary Anesthesia Machine

*The M2000 is a real workhorse, with excellent expansion, storage options...and our industry leading ten year warranty!*



*Note: Shown with optional Isoflurane Vaporizer & E-tank manifold*

The Supera M2000 comes with these features:

- Stable “A” frame design
- Large heavy duty shelf for monitors & ventilators
- Proven safety pop off valve - it saves patient lives!
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with easy change “One Seal” canister
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)
- Easily add universal control arm for Bain systems
- Built in breathing circuit and bag organizer
- Easy mounting for waste gas evacuation systems

Breathing circuits and bags included:

Adult & Pediatric Unilimb rebreathing circuit

Non-rebreathing circuit

0.5, 1, 2, 3 liter heavy duty bags

Stainless Steel IV Pole



Includes our proven safety pop off valve

# SUPERA

ANESTHESIA  
INNOVATIONS

**BROWARD A&C**  
MEDICAL GAS SPECIALISTS

954-725-1470 Ext. 403

## Features & Options

Large 22 1/2" X 13" monitor and ventilator shelf.  
Built in syringe holder and accessory mounting points



Built in circuit and bag holders help  
with moisture control.



Optional  
Waste Gas Interface  
P/N EVC630 (shown)  
*For fan evacuation systems*  
P/N EVC629 (not shown)  
*For vacuum pump systems*

Built in accessory  
mounting points

IV Pole included



Optional  
Universal Control Arm  
P/N MA2013  
*For safer non-rebreathing  
procedures*

**Includes:** Safety pop off valve,  
60" long Bain circuit and  
Supera 1 liter heavy duty  
rebreathing bag

