

## M2500 Wall Mount Veterinary Anesthesia Machine

*Full 180 degree movement with 28" super strong extension arms. Perfect for split use between multiple wet & dry tables. Great for radiology or dental areas too!*



**M2500 Includes the patented safety pop off valve**

### The Supera M2500 comes with these features:

- Articulating arm extends 28" from the wall
- Patented safety pop off valve - it saves patient lives!
- 0-4 Lpm oxygen flowmeter with mechanical stop
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with easy change "One Seal" canister
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)

### Breathing circuits and bags included:

- Adult Unilimb rebreathing circuit
- Pediatric Unilimb rebreathing circuit
- Non-rebreathing circuit
- 0.5, 1, 2, 3 liter heavy duty bags

**Note: Vaporizer & E-tank manifold are not included. Sold Separately**