

Model M3000

Non-Rebreathing Table Top

Widely used for rodents and other small animals less than 15 lbs. It takes up less than one square foot of table space and is light enough to carry between treatment areas or in the field.



Note: Vaporizer not included. Sold Separately



Optional Universal Control Arm
Manifold for Bain Circuits
P/N NA2013



The Supera M3000 comes with these features:

- 0-4 Lpm oxygen flowmeter with mechanical stop control
- Oxygen flush with safety restrictor and button shield
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)
- Balanced carrying handle
- Multiple mounting points for accessories

SUPERA

ANESTHESIA
INNOVATIONS

BROWARD & C
MEDICAL GAS SPECIALISTS

954-725-1470 Ext. 403