

# Model M4000 Rebreathing Table Top Veterinary Anesthesia Machine

Our most compact anesthesia machine. Ideal for table top, wet table or mobile carts. Also great for mobile clinics!



The M4000 comes with these features:

- Patented safety pop off valve - it saves patient lives!
- Oxygen flush with safety restrictor and button shield
- CO2 absorber with easy change "One Seal" canister
- Easy to use non-rebreathing circuit common outlet
- Syringe holder (for inflating E/T cuffs)
- Dimensions: 17.5"W x 9"L x 14"H

Note: Vaporizer not included. Sold Separately



Includes the patented safety pop off valve



I.V. Pole included



Optional Universal  
Control Arm  
Model: MA2013

# SUPERA

## ANESTHESIA INNOVATIONS

**BROWARD & C**  
MEDICAL GAS SPECIALISTS

954-725-1470 Ext. 403