



Specifications

Capacity: Up to 6 Users

Warranty: 10 Years

Flow Rate: 85 CFM

Noise Level: 61 dB

Vacuum ("Hg): 8.5" Hg

Chassis Footprint: 20" x 20" x 17"

Chassis Weight: 85 lbs

Power: 2.4 HP

Voltage: 230 Volts. (115 volts optional)

Current: 10 Amps

Tank size: 40" Tall / 12 Gal.

Tank Weight: 43 lbs

Completely Water-Less

Completely Oil-Less